





**FormaAid<sup>®</sup>** Collagen Membrane  
DOH-MD No. 003195

## GTR Membrane With Superb Stretch & Tear-resistance

### Application

- 1 Guided Tissue Regeneration (GTR)
- 2 Dental implant and oral surgery
- 3 Controlled resorption time (barrier function) and flexibility

### Feature



## Advantage

- 1 FormaAid Collagen Membrane is made of highly purified type I collagen from bovine tendon tissue.
- 2 Excellent tissue integration and barrier function for optimal bone regeneration.
- 3 Structural integrity for easy handling and suturing.
- 4 Complete resorbable, avoids possible second surgery for removing implanted membranes.
- 5 Conformable and repositionable for precise adjustment and placement.

## Specification



● 15 mm x 20 mm



● 25 mm x 30 mm




● 30 mm x 40 mm

Catalogue	Specification (1 unit/box)
FA-T1520	15 mm x 20 mm x 0.3 mm
FA-T2530	25 mm x 30 mm x 0.3 mm
FA-T3040	30 mm x 40 mm x 0.3 mm
FA-1520	15 mm x 20 mm x 0.5 mm
FA-2530	25 mm x 30 mm x 0.5 mm
FA-3040	30 mm x 40 mm x 0.5 mm

## About Us

25 years of dedication to the research and development of hyaluronic acid & collagen medical devices and health care products



**MBI** MAXIGEN BIOTECH INC.

2F, 3F, 4F, 5F, No.88, Keji 1st Rd., Guishan District, Taoyuan City 33383, Taiwan

TEL : +886-3-3287-222 FAX : +886-3-3287-333



For medical professionals reference only